RIDASCREEN[®] Histamin

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RIDASCREEN[®] Histamin

Brief information

RIDASCREEN[®] Histamin (Art. No.: R1601, 96 wells / Art. No.: R1604, 48 wells) is a competitive enzyme immunoassay for the quantitative analysis of histamine in food.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 96 or 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: different depending on the sample (see 8.)

Time requirement: sample preparation (for 10 samples)..... approx. 15 min
test implementation (incubation time)..... 90 min
test implementation (incubation time with shaker) . 70 min

Detection limit: sparkling wine, white and red wine..... 250 ppb
milk..... 100 ppb
cheese, fresh and canned fish 2.5 ppm
fish meal..... 100 ppm

Recovery rate: sparkling wine, white and red wine..... 95 %
milk..... 109 %
cheese, fresh and canned fish 100 %
fish meal..... 96 %

Specificity: Histamine 100 %
3-Methyl-histamine..... approx. 0,01 %
Tyramine n.d.
L-Phenylalanine..... n.d.
L-Histidine n.d.
L-Tyrosine n.d.
Tryptamine n.d.
5-Hydroxy-Indole-Acetic Acid n.d.
Serotonin n.d.

Used symbols:

	Contains sufficient for <n> tests		Manufacturer
	Storage temperature	LOT	Batch code
	Use by	CONT	Content
	Read entire protocol before use	REF	Art. No.
IVU	Only for in vitro use	REC	Reconstitute with
LYO	Lyophilisate		

1. Intended use

RIDASCREEN® Histamin is a competitive enzyme immunoassay for the quantitative analysis of histamine in white and red wine, sparkling wine, milk, cheese, fresh fish, canned fish and fish meal.

2. General

Histamine is a product of the decomposition of histidine caused by the growth of certain bacteria in protein rich food like fish, cheese, meat as well as sparkling wine, wine and beer. The amount of histamine formed depends on the bacterial species, the temperature and time of exposure, and may exceed 1000 ppm (mg/kg). Fish containing high levels of histamine have been associated with poisoning commonly referred to as “scombroid poisoning”. Quality control measures have to be taken to minimize the occurrence of histamine contaminated fish. Good quality fish contains less than 10 ppm histamine. The limit value for fish and products thereof is 100 mg/kg, in some countries even 50 mg/kg. The Swiss limit value for histamine in wine is 10 mg/l.

3. Test principle

After sample preparation histamine is derivatized quantitatively by an acylation reagent into N-acylhistamine. In a competitive ELISA, free acylated histamine and bound histamine compete for the antibody binding sites.

After washing, secondary antibodies labeled with peroxidase (enzyme conjugate) are added. These antibodies bind to the antibody-histamine complexes. Any unbound enzyme conjugated antibody is then removed in a washing step. Enzyme substrate and chromogen are added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue product.

The addition of the stop solution leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorption is inversely proportional to the histamine concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 96 (R1601) or 48 (R1604) measurements (including standard analyses). Each test kit contains:

Reagents for sample acylation:

Acylation plate	Acylation plate (uncoated) with 96 wells for acylation
Acylation tubes	96 x (R1601) or 48 x (R1604) Acylation tubes (plastic) for acylation
Standard 1-6	6 x Histamine standards , (4 ml each) 0 ppb, 0.5 ppb, 1.5 ppb, 5 ppb, 15 ppb, 50 ppb
Control 1	Control 1 (4 ml), concentration see vial
Control 2	Control 2 (4 ml), concentration see vial
Acylation reagent	2 x (R1601) or 1 x (R1604) Acylation reagent (1.5 ml), ready to use
Acylation buffer	Acylation buffer (22 ml), ready to use

Reagents for the enzyme immunoassay:

Microtiter plate	Microtiter plate with 96 or 48 wells (12 strips, R1601 or 6 strips, R1604 with 8 wells each) coated with histamine
Conjugate	Conjugate (12 ml, R1601 or 6 ml, R1604) peroxidase conjugated antibody, ready to use
Anti-histamine antibody	Anti-histamine antibody (12 ml, R1601 or 6 ml, R1604) ready to use
Substrate-/Chromogen	Substrate-/Chromogen (12 ml, R1601 or 6 ml, R1604) contains tetramethyl-benzidine, ready to use
Stop solution	Stop solution (12 ml), contains 0.5 N sulfuric acid, ready to use
Washing buffer	Washing buffer (20 ml), as a 50x concentrate

5. Materials required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- laboratory mincer / grinder, mortar, or homogenisator
- centrifuge and centrifugal plastic vials
- shaker, e.g. vortex mixer
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- to shorten incubation time: microtiter plate (MTP) shaker (shaking amplitude 3 mm, approx. 600 rpm)

5.2. Reagents:

- distilled or deionized water
- for milk samples: 10 mM PBS, 0.05 % Tween 20 (e.g. Sigma P3563)

6. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit reagents.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46°F).

The substrate/chromogen solution is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiration date (see  on the kit label).

Do not interchange individual reagents between kits of different lot numbers.

7. Indication of instability or deterioration of reagents

- a blue coloration of the substrate-/chromogen solution before pipetting
- an absorption smaller than 0.9 ($A_{450\text{ nm}} < 0.9$) for standard 1

8. Preparation of Samples

Acylation tubes or reaction wells of acylation plate are for single use only to avoid contaminations.

For sample preparation use only plastic vials!!!

8.1. Sparkling wine, white wine, red wine

- dilute 20 µl wine with 10 ml of distilled water
- apply 100 µl per acylation tube or per well in the acylation plate

8.2. Milk

- centrifuge 4 ml milk for 10 min / 3000 g at room temperature (20 - 25 °C / 68 - 77 °F) and remove the lipid layer
- mix 20 µl defatted milk with 4 ml PBS buffer (see 5.2.)
- apply 100 µl per acylation tube or per well in the acylation plate

8.3. Cheese, fish (fresh and canned fish)

- homogenize 10 g of sample
- add 9 ml of distilled water to 1g of homogenate and mix well
- centrifuge: 5 min / 2500 g at room temperature (20 - 25 °C / 68 - 77 °F)
- remove lipid layer or cut through lipid layer
- mix 1 ml of supernatant with 9 ml of water and mix well
- dilute 200 µl of this solution with 9.8 ml of distilled water
- apply 100 µl per acylation tube or per well in the acylation plate

8.4. Fish meal

- mix 1 g of fish meal with 200 ml of distilled water and stir for 15 min
- centrifuge aliquot: 5 min / 2500 g at room temperature (20 - 25 °C / 68 - 77 °F)
- dilute 200 µl of the supernatant with 200 ml of distilled water
- apply 100 µl per acylation tube or per well in the acylation plate

9. Test implementation

9.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **washing buffer** is provided as a 50fold concentrate. The buffer concentrate (20 - 25 °C / 68 - 77 °F) is diluted 1:50 (1+49) with distilled water before use (i.e. 10 ml buffer concentrate + 490 ml distilled water). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. four weeks.

The **acylation reagent** is provided ready to use, it has a freezing point of 18.5 °C (65 °F). When being used it must be ensured that the acylation reagent has reached room temperature and forms a homogeneous, crystal-free solution. The acylation reagent can be stored at room temperature (20 - 25 °C/ 68 - 77 °F) or brought to room temperature before use.

9.2. Test procedure for the acylation

Pipette carefully and avoid splashes.

1. All standards, controls and samples should be run in duplicate (**Art. No.: R1604, 48 wells**: Use only 48 - not all 96 - wells of the acylation plate, because the ELISA-plate contains only 48 wells). Record standard and sample positions.

Alternatively, insert required number of acylation tubes into the acylation plate and use as reaction vials.

2. Add 100 µl of each standard solution, control or prepared sample to acylation tubes or separate wells of the acylation plate.
3. Add 25 µl of the acylation reagent to each acylation tube or well.
4. Add 200 µl of the acylation buffer to each acylation tube or well, mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Use 25 µl for the ELISA.

9.3. Test procedure for ELISA

Pipette carefully and avoid splashes. Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

1. Insert a sufficient number of wells into the microwell holder for all standards, controls and samples to be run in duplicate. Record standard and sample positions.
2. Add 25 µl of acylated standard solution, control or prepared sample to separate wells.
3. Add 100 µl of the anti-histamine antibody solution to each well, mix gently by shaking the plate manually and incubate for 40 min at room temperature (20 - 25 °C / 68 - 77 °F). Alternatively, 30 min with a MTP-shaker (approx. 600 rpm)
4. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl of washing buffer (see 9.1.) and pour out the liquid again. Repeat two more times.
5. Add 100 µl of the conjugate solution to of each well. Mix gently by shaking the plate manually and incubate for 20 min at room temperature (20 - 25 °C / 68 - 77 °F). Alternatively, 10 min with a MTP-shaker (approx. 600 rpm)

6. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl of washing buffer (see 9.1.) and pour out the liquid again. Repeat two more times.
7. Add 100 µl of the substrate-/chromogen solution to each well. Mix gently by shaking the plate manually and incubate for 15 min at room temperature (20 - 25 °C / 68 - 77 °F) **in the dark**. Alternatively, 15 min with a MTP-shaker (approx. 600 rpm).
8. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm against an air blank. Read within 10 min after addition of stop solution.

10. Results

A special software, the RIDA[®]SOFT Win (Art. No. Z9999), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays.

The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

The histamine concentration in µg/kg (ppb) corresponding to the absorbance of each sample is read from the calibration curve and then further multiplied by the corresponding dilution factor. The dilution factor depends on the sample preparation and is listed in the table below.

Sample preparation	Sample	Dilution factor
8.1.	red wine, white wine and sparkling wine	500
8.2.	milk	200
8.3.	cheese, fresh and canned fish	5 000
8.4.	fish meal	200 000

Recommendation

In order to ensure a high analytical performance:

- further dilute samples with absorbances < 0.4 and analyse again
- analyze each sample in duplicate
- the controls should be measured with a maximum deviation of +/- 30 % (otherwise repeat the assay)
- the optimal temperature for using the assay is 20 - 25 °C / 68 - 77 °F.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG