RIDA[®] QUICK Histamin

Colorimetric assay for the quantitative determination of
histamine in fresh fish and fish meal

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www.r-biopharm.com

For further questions please contact your local distributor.

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RIDA[®] and RIDASCREEN[®]
are registered trademarks of R-Biopharm AG
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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Brief information

RIDA[®]QUICK Histamin (Art. No.: R1603) is a colorimetric assay for the quantitative determination of histamine in fresh fish and fish meal.

The test kit is sufficient for 48 determinations (sample determinations).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation: Extraction, clean up with an ion-exchange column

Time requirement: Sample preparation (for 10 samples) approx. 15 min
 Test implementation (incubation time)..... approx. 10 min

Detection limit: Fresh fish20 ppm
 (corresponding to the Fish meal125 ppm
 standard substance)

Cross-reactions: The specificity of the RIDA[®]QUICK Histamin was determined by analyzing the cross-reactivities to corresponding substances in buffer system.

Histamin	100 %
L-Histidine	< 1
Tyramine	< 1
L-Tyrosine	< 1
L-Phenylalanine.....	n.d.
L-Tryptophan	n.d.
Non-aromatic amino acids.....	n.d.
Tryptamine	n.d.
Serotonine	n.d.

Limitation: Fish such as tuna contain naturally very high levels of histidine concentrations which disturb the assay.

Used symbols:

	Contains sufficient for <n> tests		Manufacturer
	Storage temperature	LOT	Batch code
	Use by	CONT	Content
	Read entire protocol before use	REF	Art. No.
IVU	Only for in vitro use		

1. Intended use

RIDA[®]QUICK Histamin is a colorimetric assay for the quantitative analysis of histamine in fish meal. If fresh fish samples are to be analyzed a validation has to be carried out (employ histamine-free and histamine containing (spiked) samples for each fish species).

2. General

Histamine is one of the biogenic amines. The amount of histamine formed depends on the bacterial species, the temperature and time of exposure, and may exceed 1000 ppm (mg/kg). If concentrations of amines exceed normal levels, possibly due to bacterial contamination of food, harmful effects may occur. The contents of histamine in food is regarded as criterion of quality.

3. Test principle

The RIDA[®]QUICK Histamin assay is a colorimetric assay. After extraction of the sample, a clean-up is carried out with an ion exchange column to separate histamine from impurities. After separation histamine is derivatized and the resulting dye is determined in a photometer at 450 nm wavelength. The optical density of the dye is proportional to the concentration of histamine in the sample. A standard curve allows the quantification of the histamine concentration in the sample.

4. Reagents provided

Each kit contains sufficient materials for 48 samples. Each test kit contains:

6 x 8	Ion exchanger	Ion exchanger columns
1 x	Microtiter plate	Microtiter plate 96 wells non-coated (12 strips with 8 removable wells each)
6 x	Standard N	Histamine standard solutions, 0.9 ml each *) 0 ppm (zero standard), 2.5 ppm, 5 ppm, 10 ppm, 20 ppm, 40 ppm, ready to use
1 x	Elution buffer	Elution buffer (50 ml), ready to use
1 x	Washing buffer	Washing buffer (50 ml), 10 fold concentrate
1 x	Color reagent 1	Color reagent 1 (5.4 ml)
1 x	Color reagent 2	Color reagent 2 (5.4 ml)

*) The dilution factor 4 of the standards during the color reaction has already been taken into account.

5. Materials required but not provided

5.1. Equipment

- microtiter plate or strip photometer (450 nm)
- grinder (mill)
- graduated cylinder: 100 ml plastic or glass
- plastic tubes for preparing sample extract
- optional: shaker
- vortex
- filter paper: Whatman No. 1 or equivalent
- variable 20 - 200 µl and 200 - 1000 µl micropipettes

5.2. Reagents

- isopropanol
- distilled or deionized water
- 70 % isopropanol solution: prepare 70 % of isopropanol solution by mixing 70 ml pure isopropanol with 30 ml of distilled or deionized water

6. Warnings and precautions for the users

This test should be performed only by trained laboratory staff. The instructions for performing the test must be followed strictly.

Isopropanol is inflammable.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused ion exchanger columns to their original foil bag and reseal them.

No quality guarantee is accepted after the expiration date (see  on the kit label).

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

Precipitates in any of the reagents may indicate spoilage and reagents may be deteriorated. Replace test kit if precipitates occur.

9. Preparation of Samples

9.1. Preliminary comments

1. Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.
2. Return all reagents to 2 - 8 °C (36 - 46 °F) immediately after use.
3. The washing buffer is provided as a 10fold concentrate. Before use dilute the washing buffer 1:10 (1+9) with distilled water, e. g. 1 ml buffer with 9 ml distilled water.
4. Notes regarding clean up procedure: Displace the top stopper of the column first, then remove the bottom stopper to avoid splashing off the gel matrix. The sample as well as the washing buffer should be applied carefully on the top of the column in order to avoid flushing the gel material.
5. Do not allow the columns to dry between working steps.
6. Reproducibility is largely dependent upon the consistency with which the columns are washed. Carefully follow the recommended washing sequence as outlined.

9.2. Fresh fish

Extraction

- homogenize 50 g of fish with 50 ml of distilled water
- transfer 1 g of homogenate into a polypropylene tube
- add 3 ml of isopropanol (100 %)
- mix with a vortex
- centrifuge: 3 min / at 2000 g / room temperature (20 - 25 °C / 68 - 77 °F)
- use 1 ml of the supernatant for the clean up procedure

Clean up

- precondition the ion exchange column with 3 ml of diluted washing buffer (see 9.1.)
- add 1 ml of the supernatant onto the column and allow it to flow slowly into the column
- wash the column with 1 ml of 70 % isopropanol solution
- wash the column twice (2 x) with 3 ml of diluted washing buffer (see 9.1.)
- discard the washing solution and use a new tube to collect the histamine eluate
- pipet 0.5 ml of elution buffer carefully onto the column, allow it to flow slowly into the column and pipet again 0.5 ml of elution buffer onto the ion-exchanger column
- press some ml of air through the column to obtain the entire eluate
- mix the eluate with a vortex
- use 200 µl per well of the mixed eluate for the color reaction

9.3. Fish meal

Extraction

- weigh 2 g of ground sample and add it to a suitable container with 20 ml of 70 % isopropanol solution
- shake vigorously for four minutes (manually or with shaker)
- filter the extract through a paper filter (e. g. Whatman No. 1 filter), centrifuge or just let the solids settle for a few minutes
- use 200 µl of the filtrate or supernatant for the clean up procedure

Clean up

- precondition the ion exchange column with 3 ml of diluted washing buffer (see 9.1.)
- add 200 µl of the supernatant onto the column
- wash the column twice (2 x) with 3 ml of diluted washing buffer (see 9.1.)
- discard the washing solution and use a new tube to collect the histamine eluate
- pipet 0.5 ml of elution buffer carefully onto the column, allow it to flow slowly into the column and pipet again 0.5 ml of elution buffer onto the ion-exchanger column
- press some ml of air through the column to obtain the entire eluate
- mix the eluate with a vortex
- use 200 µl per well of the mixed eluate for the color reaction

10. Test procedure for the colorimetric assay

1. Mix the color reagents 1 and 2 in equal amounts (1+1). This mixture is stable for at least 10 hours.
2. Insert a sufficient number of wells into the microwell holder for all standards and samples to be run. Record standard and sample positions.
3. Pipet 50 µl of standard and 150 µl elution buffer to each standard well. Use a new pipette tip for each standard.
4. Add 200 µl of each sample eluate to a new well.
5. Add 50 µl of the color reagent mixture to each well (samples and standards)
6. Measure the absorbance at 450 nm against an air blank. Read within 10 min.

11. Results

A special software, the RIDA®SOFT Win (Art. No. Z9999), is available to evaluate the RIDA®QUICK Histamin assay.

In order to obtain the histamine concentration in mg/kg (ppm) actually contained in a sample, the concentration read from the calibration curve must be further multiplied by the corresponding dilution factor. When working in accordance with the described procedure (see 9.), the dilution factor is as follows:

fresh fish.....	8
fish meal.....	50

Therefore, the histamine measuring range of the standard curve is between 20 and 320 mg/kg (ppm) for fresh fish and between 125 to 2000 mg/kg (ppm) for fish meal.

Limitations:

- Very high concentrations of histidine and biogenic amines disturb the assay.

Important Note:

- Fresh fish samples have to be validated: Employ histamine-free and histamine containing (spiked) samples for each fish species. A spiking procedure and a spiking solution is available from your local distributor.
- Extract each sample twice.
- Further dilute samples with an absorbance > 3 and analyze again. The same sample extract can be used for retesting in applying less sample extract onto the column. If e.g. only 100 µl are applied, the dilution factor is 100 (the corresponding measuring range of the standard curve is 250 to 4000 mg/kg).
- This test should only be carried out by trained laboratory staff.

R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, R-Biopharm will provide a replacement product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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