RIDASCREEN[®] Histamine (enzymatic)

Art. No. R1605

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RIDASCREEN® Histamine (enzymatic)

Brief information

RIDASCREEN® Histamine (enzymatic) (Art. No. 1605) is an enzymatic test in microtiter plate format for the quantitative determination of histamine in fresh fish, canned fish, fish meal, wine (in combination with Art. No. R1699), cheese and milk.

All reagents required for the enzyme assay - including standards - are contained in the test kit.

The test kit is sufficient for 96 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization and extraction
Time requirement:	sample preparation (10 fish samples. approx. 25 min test implementation (incubation time) 15 min
Standard:	histamine
Limit of detection:	fresh fish, canned fish 0.75 mg/kg fish meal 3.75 mg/kg wine 0.54 mg/kg cheese 0.75 mg/kg
Limit of quantification:	fresh fish, canned fish 2 mg/kg fish meal 10 mg/kg wine 1.4 mg/kg cheese 2 mg/kg
Recovery:	approx. 85 - 105 %
Specificity:	The enzyme is specific for histamine. There are no side reactions to L-histidine, L-tryptophan, serotonin, L-tyrosine, tyramine, L-lysine, trimethylamine.

The specificity of RIDASCREEN® Histamine (enzymatic) has been evaluated by the determination of the side activities of the respective substances in a buffer system (see validation report).

Related products

RIDA® Sample Decolorant (R1699) (for wine analysis)

1. Intended use

RIDASCREEN® Histamine (enzymatic) is used for the quantitative determination of histamine in fresh fish, canned fish, fish meal, wine, cheese and milk.

2. General

Histamine is a product of the decomposition of histidine caused by the growth of certain bacteria in protein rich food like fish, cheese and wine.

The amount of histamine formed depends on the bacterial species, the temperature and the time of exposure, and may exceed in fish 1000 mg/kg. Fish containing high levels of histamine have been associated with poisoning commonly referred to as “scombroid poisoning”.

Fish of good quality contains less than 10 mg/kg histamine. The limit value for fish and products thereof is between 50 and 200 mg/kg depending on the fish species and the country. Similar high concentrations can be found in cheese, while histamine in wine usually does not exceed 15 mg/l.

3. Test principle

The basis is an enzymatic reaction. The microtiterplate is coated with a reagent (electron carrier) and a dye.

Histamine-dehydrogenase catalyzes the oxidation of histamine to imidacetaldehyde in presence of an electron carrier and a dye. The forming dye is measured at 450 nm and is proportional to the histamine concentration.

4. Reagents provided

Each kit contains sufficient materials for 96 measurements (including standard analyses). All reagents are ready to use. Each test kit contains:

Component	Cap color	Format		Volume
Microtiter plate	-	Ready to use		96 wells
Buffer	Transparent	Ready to use		15 ml
Standard 1	White	Ready to use	0 mg/l	1.2 ml
Standard 2	White	Ready to use	1 mg/l	1.2 ml
Standard 3	White	Ready to use	5 mg/l	1.2 ml
Standard 4	White	Ready to use	10 mg/l	1.2 ml
Standard 5	White	Ready to use	15 mg/l	1.2 ml
Standard 6	White	Ready to use	20 mg/l	1.2 ml
Enzyme solution	Red	Ready to use		1 ml
Spiking solution	blue	Ready to use	10 mg/ml	1 ml

5. Reagents required but not provided

5.1. Equipment:

For the sample preparation only plastic vials (no glass) should be used.

- vortex
- variable 100 - 200 µl micropipettes
- multistepper-pipette with 10 µl, 100 µl, 150 µl, 200 µl pipetting volume (e.g. Eppendorf combitips)
- microtiter plate spectrophotometer (450 nm)
- RIDA®SOFT Win.net (Z9996); use **version 1.92** or newer for calculation

Fish analysis

- 50 ml PP screw-cap vials (e.g. Greiner; Product No. 227261)
- homogenizer: laboratory mincer / grinder, Ultra-Turrax or mixer
- water bath (100 °C / 212 °F)
- centrifuge, centrifuge vials

Wine and cheese analysis

- 1.5 - 2 ml plastic vials (e.g. Eppendorf tube)
- mixer, blender (e.g. Mr. Magic)
- microcentrifuge

5.2. Reagents:

- distilled or deionized water

Wine analysis

- RIDA® Sample Decolorant (R1699) for sample preparation

Cheese analysis

- 1 N perchloric acid (e.g. Bernd Kraft, Art. No. 20126.2700)
- 1 N KOH (e.g. Carl Roth, Art. No. K017.1)

6. Warnings and precautions for the users

This test should be carried out only by trained laboratory technicians. The instructions for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

No quality guarantee is accepted after the expiration date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

If the difference between absorption of $A_2 - A_1$ of standard 6 is differing strongly from the certificate of analysis, then a deterioration of reagents is likely.

9. Preparation of samples

The histamine concentration increases also in refrigerated food, therefore it is important to carry out the laboratory analysis directly after sampling. **Use only plastic vials for the sample preparation.**

9.1. Preliminary comments

Dirty laboratory equipment can lead to contamination of samples and test components with histamine.

9.2. Fresh fish and canned fish

- weigh 5 g of homogenized sample into a 50 ml PP-screw cap vial
- pour 20 ml of dist. water over the sample.
- close the vial and shake or vortex until the sample is evenly suspended
- incubate the sample for 20 min in a boiling water bath and shake from time to time
- allow the sample to cool
- centrifuge for 10 min at least 2500 g (if possible at 4 °C (39 °F)
Alternatively transfer 2 ml of the extract into centrifugation vials and centrifuge at high speed (> 10000g) for 2 min in a micro-centrifuge.
If a layer of fat has formed on top of the extract after centrifugation then carefully penetrate through the fat layer and pipette the sample into a new vial and centrifuge again. (If the extract is cloudy then filter with a syringe filter.
- pipette 100 µl of the undiluted **clear extract** in the assay).

9.3. Fish meal

- weigh 1 g of homogenized fish meal into a 50 ml PP-screw cap vial
- pour 25 ml of dist. water over the sample
- close the vial and shake or vortex until the sample is evenly suspended
- incubate the sample for 20 min in a boiling water bath and shake from time to time
- allow the sample to cool
- centrifuge for 10 min at least 2500 g (if possible at 4 °C (39 °F)
Alternatively transfer 2 ml of the extract into centrifugation vials and centrifuge at high speed (> 10000g) for 2 min in a micro-centrifuge.
If a layer of fat has formed on top of the extract after centrifugation then carefully penetrate through the fat layer and pipette the sample in a new vial and centrifuge again. (If the extract is cloudy then filter with a syringe filter).
- Pipette 100 µl the undiluted **clear extract** in the assay.

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days or at -20 °C (-4 °F) for three months.

9.4. Wine

Wine (red wine, white wine, Rosé) contains anthocyanins and polyphenols which have to be removed prior to the enzymatic analysis.

- Add 200 µl Reagent 1 of RIDA[®] Sample Decolorant to 200 µl wine (rinse pipette tip)
- Add 200 µl Reagent 2 of RIDA[®] Sample Decolorant, vortex and incubate 5 min at room temperature, then centrifuge 2 min at a minimum of 14000 g (room temperature)
- Transfer 500 µl of the supernatant in a new vial
- Add 100 µl Reagent 3 of RIDA[®] Sample Decolorant, vortex and incubate 5 min at room temperature then centrifuge at 14000 g (room temperature)
- Use 100 µl of the **clear extract** in the assay

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 1 day.

9.5. Cheese

If the cheese rind is edible then it should also be part of the representative sample.

- Weigh in 60 g of cheese into the plastic cup of blender and add 140 ml of water. Mix for 15 sec. to suspend all particles. Repeat two more times.
- Remove the fat layer with a spatula. Transfer 50 ml of the unclear solution to an plastic vial.
- Centrifuge at 3500 g for 10 min (at 4°C if possible).
- Remove the fat layer with a spatula and transfer 0,8 ml of the clear solution to an plastic vial.
- Add 200 µl 1 N perchloric acid, mix well.
- Centrifuge at 20000 g for 2 min at room temperature.
- Transfer 500 µl of the clear supernatant to a new Eppendorf centrifuge vial and add 100 µl 1 N KOH
- Mix and centrifuge 2 min at least 20000 g at room temperature
- Use 100 µl of **clear extract** per

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 1 day.

10. Test implementation

10.1. Test preparation

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) and mix carefully before use.

10.2. Test implementation

Insert a sufficient number of wells into the microwell holder for all standard and samples to be run in duplicate. Record standard and sample positions.

1. Using a multi-stepper, add 150 µl of buffer to the wells and **carefully shake the plate manually**.
2. Add 100 µl of standards, controls or prepared samples to separate wells in duplicate. Flush the pipette tips. **Thereafter, carefully shake the plate manually**.
3. After 3 min the absorption A1 is measured at 450 nm without a reference filter.*
4. Add 10 µl of the blue dyed enzyme solution with a multi-stepper to each well. Thereafter, **carefully shake the plate manually**.
5. After 10 min the absorption A2 is measured at 450 nm without a reference filter. For evaluation use the RIDA®SOFT Win.net, select A1 and A2 and calculate.

*For **fish** and **wine** a creep reaction may occur during the subsequent enzymatic analysis: If during the 3 minute incubation an **increasing yellow coloration** is observed, then test for ascorbic acid concentrations lower than 250 mg/l. If no ascorbic acid is observed then increase the incubation time to 10 min.

11. Results

At R-Biopharm AG a special software, the RIDA[®]SOFT Win.net (Art. Nr. Z9996, use version 1.88 or newer), is available for evaluation of the RIDASCREEN[®] assays. The calculation should be done using linear regression, the dilution factor of the samples has to be taken into account. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

Sample	Dilution factor	Measuring range
Fresh/canned fish	5	2 - 100 mg/kg
Fish meal	25	10 - 500 mg/kg
Wine	3.6	1.4 - 72 mg/kg
Cheese	5	2 - 100 mg/kg

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of histamine contamination.

A further dilution with water and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 6.

Recommendation:

In order to ensure a high analytical performance:

- Each sample should be measured in duplicate.
- Employ also histamine-free and histamine-containing (spiked) samples as test controls.

To verify the correct and interference-free implementation of the test procedure the spike solution contained in the kit should be used for spiking.

10 mg/l spike control

- 1:1000 dilution of the Spiking Solution Histamin (10 mg/ml) in two steps e.g.
100 µl + 3.9 ml water (1:40)
200 µl thereof + 4.8 ml water (1:1000 in total)
- Use 100 µl in the assay.

50 mg/kg histamine spiked fish sample

- pipette 25 µl of the Spiking Solution Histamine (10 mg/ml) to 5 g homogenized fish
- extract the sample with 20 ml of water (see 9.2.)

250 mg/kg histamine spike fish meal sample

- pipette 25 µl of the Spiking Solution Histamine (10 mg/ml) to 1 g homogenized fish meal sample
- extract the sample with 25 ml of water (see 9.3.)

10 mg/l histamine spiked wine sample

- pipette 10 µl of the Spiking Solution Histamine (10 mg/ml) to 10 ml wine
- extract the sample with RIDA[®] Sample Decolorant (see 9.4)

Please note:

For further product information and applications (e.g. milk) please contact your local distributor or R-Biopharm AG.

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG