

RIDASCREEN® FAST Peanut

Art. No. R6202

Product information





RIDASCREEN® FAST Peanut

Art. No. R6202

General

Although peanut (*Arachis hypogaea*) belongs to the legumes and not to the nuts, its taste and its usual consumption is similar to tree nuts. So to speak, peanut is one of the most popular “nuts” worldwide. Peanuts are very common ingredients of chocolate products and are also consumed as snacks and in bakery products such as cookies.

Peanut is also one of the five foods which are the most frequent cause of food allergies. Allergies against peanuts are among the most severe allergies and can lead to life-threatening anaphylactic shocks.

Consumption of probably less than a milligram of peanut can induce allergic reactions in very sensitized individuals. Therefore, labelling of peanut containing products is very important and compulsory in the USA and the European Union as well as in many other countries.

The main threat for sensitized individuals is contamination of peanut in food which is intended to be peanut-free (e.g. peanut traces in chocolate which should not contain peanuts). Therefore, it is important to test foods which may have been contaminated during their production process for the presence of unintended traces of peanuts.

Sample preparation

Extraction

The sample extraction is carried out as described in the instruction for use. For tannin and polyphenol containing samples such as spices and chocolate, the addition of skim milk powder to the sample is recommended.

Additionally, the addition of skim milk powder is also recommended for the extraction of ice cream and caramel, as described in the instruction for use.

Spiking instruction

Detailed information on spiking of samples can be found in appendix 2.

Calibration curve

The basis of the sandwich ELISA is an antigen-antibody reaction. Detailed information can be found in the instructions for use. As a calibrator, NIST peanut butter 2387 with a protein content of 22.2 % is used. The result is expressed as mg/kg peanut.

A typical standard curve is shown in appendix 1. To generate the calibration curve, a cubic spline function is recommended.

Specificity

The antibodies specifically detect antigens from peanuts, mainly allergenic proteins such as Ara h 1 to Ara h 3. Cross-reactivities to green peas, lentils and fenugreek were observed.

To recheck the cross-reactivity we recommend using a non-immunoassay method such as PCR, SureFood® ALLERGEN ID.

To generate the data presented in this validation report, equipment as described in appendix 3 was used.



RIDASCREEN® FAST Peanut

Art. No. R6202

Sensitivity

Limit of Detection (LOD)

Different zero matrices have been extracted 10 times each to estimate the **Limit of Detection**. A second order polynomial function was used for the extrapolation of the standard curve to calculate from the OD values of the zero matrices to the respective allergen concentration. The **Limit of Detection** was then estimated by the mean of these concentrations plus 3 times standard deviation (dilution factor for sample preparation is included).

The **Limit of Detection** is matrix-dependent.

Table 1: Estimation of the Limit of Detection (LOD).

Matrix	Standard 1 (OD*)	Standard 2 (OD*)	Mean of samples (OD*)	Mean (mg/kg)	Standard deviation mean (mg/kg)	Mean + 3x standard deviation (mg/kg)
Chocolate	0.212	0.703	0.220	0.04	0.05	0.19
Cookies	0.212	0.703	0.198	-0.07	0.08	0.16
Ice cream	0.212	0.703	0.186	-0.12	0.05	0.04
Mean value (mg/kg)						0.13

* OD = optical density

Limit of Quantification (LOQ)

The **Limit of Quantification** or the lowest concentration that can be determined in a sample with acceptable precision (repeatability) and accuracy under the stated conditions of the test was proofed to be 2.5 mg/kg (dilution factor included). The LOQ is proofed by showing the range of samples spiked at the level of standard 2. Each sample was extracted 10 times and measured in duplicates. The measured sample ranges are similar to the range of standard 2, measured 10 times. All measured coefficients of variation are within an acceptable range showing that samples with a concentration of 2.5 mg/kg can be measured with good precision and accuracy.

Table 2: Determination of the Limit of Quantification by extracting spiked samples ten times.

Sample	Measured concentration (mg/kg)	Standard deviation (SD)	Coefficient of variation (%)
Standard 2 (n=10)	2.38	0.24	9.9
Chocolate (n=10)	2.24	0.11	5.1
Cookie (n=10)	2.98	0.16	5.3



RIDASCREEN® FAST Peanut

Art. No. R6202

Variance of the ELISA

To determine the reproducibility and repeatability of the ELISA, intra-assay and inter-assay coefficients were calculated.

Intra-Assay Variation (repeatability)

Within run variation was calculated by measuring the standards in one assay run (n=6). The within assay variation over the standard curve was 8.1 % (OD-value based, data not shown). The Intra-Assay Variation of spiked samples has been found in a good range.

Table 3: Determination of the Intra-assay variation (repeatability) of spiked samples.

Matrix	Expected concentration (mg/kg)	Concentration found (mg/kg)	Coefficient of variation (%)	Recovery (%)
Cookie	5.0	4.9	5.2	98
Chocolate	20.0	21.1	10.9	104
Ice cream	0	<LOQ		
Mean value			8.0	101

Inter-Assay Variation (reproducibility)

Between run variation was determined by repeated measurements of the standard curves from different test kits of one production batch, performed by three technicians (n=3 runs). In each run duplicates of samples and standards have been measured. The mean between assay variance was 7.5 % over the standard curve (OD-value based, data not shown). The Inter-Assay Variance of spiked food samples has been found in a good range.

Table 4: Determination of the Inter-Assay Variation (reproducibility) of spiked samples.

Matrix	Expected concentration (mg/kg)	Concentration found (mg/kg)	Coefficient of variation (%)	Recovery (%)
Cookie	5.0	5.4	1.6	108
Chocolate	200	20.2	1.6	101
Ice cream	0	< LOQ		
Mean value			1.6	104



RIDASCREEN® FAST Peanut

Art. No. R6202

Cross reactivity

Food samples from different origins were evaluated as a pure commodity (100 % level) for their cross reactivity in the ELISA. Samples found well below the LOQ showed no cross reactivity. Other samples indicated in bold letters showed a signal over LOQ, for these samples a cross reactivity was demonstrated.

Table 5: Overview of different foods tested for cross reactivity at a hundred percent level. All samples have been extracted according to the instruction for use.

Sample description	OD Standard 1*	OD Standard 2*	OD Sample	Result (mg/kg)
Cereals/pseudocereals				
Amaranth	0.206	0.659	0.149	<LOQ
Barley	0.206	0.659	0.185	<LOQ
Buckwheat flour	0.206	0.659	0.381	<LOQ
Corn flour	0.206	0.659	0.106	<LOQ
Kamut	0.206	0.659	0.140	<LOQ
Millet flour	0.206	0.659	0.087	<LOQ
Oat flour	0.206	0.659	0.091	<LOQ
Quinoa	0.206	0.659	0.076	<LOQ
Rice flour	0.206	0.659	0.075	<LOQ
Rye	0.23	0.637	0.307	<LOQ
Spelt	0.23	0.637	0.279	<LOQ
Teff flour	0.228	0.572	0.169	<LOQ
Wheat	0.228	0.572	0.232	<LOQ
Beans/lentils/peas				
Chick peas	0.228	0.572	0.268	<LOQ
Green bean	0.244	0.696	0.299	<LOQ
Green peas	0.187	0.621	0.840	3.65
Kidney beans	0.228	0.572	0.276	<LOQ
Lentils	0.187	0.621	0.948	4.35
Lima beans	0.228	0.572	0.174	<LOQ
Pinto beans	0.228	0.572	0.363	<LOQ
Soy milk	0.228	0.572	0.214	<LOQ
Soya flour	0.230	0.637	0.635	<LOQ
White beans	0.134	0.484	0.143	<LOQ
Soya flour	0.258	0.337	0.214	<LOQ
White beans	0.258	0.337	0.212	<LOQ
Seeds/stones/kernel				
Linseed	0.228	0.572	0.206	<LOQ
Pine kernel	0.228	0.572	0.122	<LOQ
Poppy seeds	0.235	0.718	0.233	<LOQ
Pumpkin kernels	0.228	0.572	0.287	<LOQ
Rape seed	0.235	0.718	0.223	<LOQ
Sesame	0.235	0.718	0.232	<LOQ
Sunflower kernels	0.235	0.718	0.174	<LOQ

Product information

RIDASCREEN® FAST Peanut

Art. No. R6202

Sample description	OD Standard 1*	OD Standard 2*	OD Sample	Result (mg/kg)
Nuts				
Brazil nut	0.235	0.718	0.170	<LOQ
Cashew	0.235	0.718	0.341	<LOQ
Cashew (roasted)	0.235	0.718	0.208	<LOQ
Chestnut	0.235	0.718	0.180	<LOQ
Macadamia	0.235	0.718	0.267	<LOQ
Pecan nut	0.235	0.718	0.182	<LOQ
Pistachio nut	0.235	0.718	0.199	<LOQ
Walnut	0.235	0.718	0.211	<LOQ
Walnut (roasted)	0.235	0.718	0.204	<LOQ
Spices				
Aniseed	0.235	0.718	0.266	<LOQ
Basil	0.214	0.590	0.493	<LOQ
Caraway	0.214	0.590	0.195	<LOQ
Cardamon	0.214	0.590	0.152	<LOQ
Celery (powder)	0.214	0.590	0.187	<LOQ
Chilli	0.214	0.590	0.211	<LOQ
Cinnamon	0.214	0.590	0.213	<LOQ
Cloves	0.198	0.682	0.320	<LOQ
Cumin	0.214	0.590	0.194	<LOQ
Curcuma	0.198	0.682	0.168	<LOQ
Curry	0.214	0.590	0.271	<LOQ
Fennel seed	0.214	0.590	0.185	<LOQ
Fenugreek	0.198	0.682	0.779	2.88
Garlic	0.214	0.590	0.531	<LOQ
Ginger	0.214	0.590	0.208	<LOQ
Marjoram	0.141	0.544	0.166	<LOQ
Mustard seeds (black-brown)	0.214	0.590	0.173	<LOQ
Mustard seeds (yellow)	0.141	0.544	0.147	<LOQ
Nutmeg	0.141	0.544	0.112	<LOQ
Onion (powder)	0.141	0.544	0.120	<LOQ
Paprika (sweet)	0.141	0.544	0.382	<LOQ
Pepper (black)	0.141	0.544	0.128	<LOQ
Pepper (white)	0.141	0.544	0.099	<LOQ
Meat				
Beef (cooked)	0.206	0.659	0.177	<LOQ
Chicken (cooked)	0.206	0.659	0.105	<LOQ
Pork (cooked)	0.206	0.659	0.085	<LOQ
Turkey (cooked)	0.206	0.659	0.125	<LOQ
Various				
Apple juice	0.141	0.544	0.118	<LOQ
Bovine gelatin	0.141	0.544	0.141	<LOQ
Bovine serum albumin	0.141	0.544	0.133	<LOQ
Cocoa (powder)	0.141	0.544	0.128	<LOQ
Coconut	0.141	0.544	0.131	<LOQ
Fetal calf serum	0.134	0.484	0.088	<LOQ
Fish gelatin	0.141	0.544	0.150	<LOQ
Orange juice	0.141	0.544	0.121	<LOQ
Porcine gelatin	0.141	0.544	0.173	<LOQ
Skim milk (powder)	0.141	0.544	0.140	<LOQ
Soya lecithin	0.141	0.544	0.220	<LOQ
Wheat starch	0.141	0.544	0.169	<LOQ
Whole egg (powder)	0.141	0.544	0.150	<LOQ

* Results are derived from different runs, i.e. slightly different absorbance values of standards were obtained

Product information

RIDASCREEN® FAST Peanut

Art. No. R6202

Due to the high potential of cross reactivity to nuts, the reactivity of the ELISA to almonds and hazelnuts was examined in detail (table 6).

Table 6: Determination of cross reactivity to roasted almonds and hazelnuts.


Sample description	OD Std 1	OD Std 2	OD Sample	Result (mg/kg)
Almond				
Untreated	0.185	0.680	0.131	<LOQ
Roasted 10 min 130 °C (266 °F)	0.185	0.680	0.126	<LOQ
Roasted 10 min 150 °C (302 °F)	0.185	0.680	0.147	<LOQ
Roasted 10 min 180 °C (356 °F)	0.185	0.680	0.129	<LOQ
Roasted 10 min 200 °C (392 °F)	0.185	0.680	0.123	<LOQ
Hazelnut				
Untreated	0.185	0.680	0.140	<LOQ
Roasted 10 min 130 °C (266 °F)	0.185	0.680	0.125	<LOQ
Roasted 10 min 150 °C (302 °F)	0.185	0.680	0.158	<LOQ
Roasted 10 min 180 °C (356 °F)	0.185	0.680	0.173	<LOQ
Roasted 10 min 200 °C (392 °F)	0.185	0.680	0.144	<LOQ

Recovery experiments

To measure the accuracy, zero samples have been spiked with an extract as described in appendix 2.

Table 7: Different spike levels for various food samples.

Declaration	Final dilution	Concentration calculated (mg/kg)	Concentration measured (mg/kg)	Recovery (%)
Chocolate, peanut free	20	0	< LOQ	
Peanut oil, peanut free	20	0	< LOQ	
Butter cookie, peanut free	20	0	< LOQ	
Fapas rounds peanut in chocolate				
FAPAS 2741 B	20	14.5	16.5	113.8
FAPAS 2752 B	20	20.2	23.4	115.8
FAPAS 2766 B	20	22.2	19.7	88.7
FAPAS 2776 B	40	40.1	39.3	98.0
Pilot plant produced cookies with defined amounts of peanut				
Cookies	20	0	< LOQ	
Cookies	20	5	5.8	116.0
Cookies	20	10	11.2	112.0
Cookies	40	20	21.7	108.5
Pilot plant produced chocolate with defined amounts of peanut (Par. 35 collab study samples)				
Chocolate, peanut free	20	0	< LOQ	
Dark chocolate	20	10	12.9	129.0
Dark chocolate	80	40	48.7	121.8
Dark chocolate	160	80	79.5	99.3
Milk chocolate	20	5	5.9	118.0
Milk chocolate	20	10	11.7	116.6
Milk chocolate	80	20	21.3	106.4
Milk chocolate	20	40	46.6	116.4
Candy	20	10	13.2	132.2



Product information

RIDASCREEN® FAST Peanut

Art. No. R6202

Table 8: Different spike levels for various food samples.

Matrix	Spike (mg/kg)	Concentration found (mg/kg)	Recovery (%)
Cereals	4.5	5.4	120.7
	5.0	5.3	105.3
	6.0	6.2	102.9
	11.1	10.9	97.7
Chocolate	5.0	4.4	88.0
	5.4	4.6	85.2
	5.5	5.4	98.2
	6.5	6.2	95.8
	10.0	9.5	94.5
Cookies	4.9	4.8	98.0
	5.0	4.6	92.0
	6.2	5.6	89.8
Ice cream	5.8	6.2	107.3
	6.0	6.1	99.2
	9.1	9.0	98.9
	10.3	9.8	95.1
	15.5	13.0	83.5
Overall mean			97.2

Spiked processed matrices

Food processing such as roasting can alter protein structures and therefore influence antibody binding which, as a consequence, can lead to reduced recoveries.

For untreated peanuts, an overestimation was measured of approx. 350 %. It should be noted that untreated peanuts are very rarely used in commercial foods. The overestimation of untreated peanuts is due to the standardization to the NIST peanut butter 2387, which contains roasted peanuts. The degree of roasting of the NIST peanut butter is also the range of roasting used for of the vast majority of peanuts in commercial foods. Roasting at higher temperatures reduces the recovery significantly. Please note that these high levels of roasting are rarely used as the peanut color changes to dark brown and the taste strongly deteriorates.

Table 9: Recovery (%) of peanuts at different degrees of roasting.

Sample	Peanut concentration (mg/kg)	Measured concentration (mg/kg peanut)	Recovery (%)
Peanut untreated	10	34.7	347
Peanut, roasted 10 min at 130 °C (266 °F)	10	10.5	105
Peanut, roasted 10 min at 150 °C (302 °F)	10	10.9	109
Peanut, roasted 10 min at 180 °C (356 °F)	10	< LOQ	< 25 (measured below LOQ)
Peanut, roasted 10 min at 200 °C (392 °F)	10	< LOQ	< 25 (measured below LOQ)

RIDASCREEN® FAST Peanut

Art. No. R6202

Spiked spices

Spices contain large amounts of polyphenols, inhibitors and disturbing substances which might have an impact on any analytical method. Therefore, a wide panel of spices as well as cocoa was spiked according to appendix 2 and analyzed for recovery.

Table 10: Recovery (%) of spiked spices and cocoa using extraction with skim milk powder.

Sample	Expected concentration (mg/kg)	Concentration found (mg/kg)	Recovery (%)
Anise seed	10	9.5	95
Basil	10	< LOQ	< 25
Caraway	10	6.0	60
Cinnamon	10	7.6	76
Cloves	10	< LOQ	< 25
Cocoa	10	5.6	56
Cumin	10	8.3	83
Curcuma	10	8.7	87
Fennel seed	10	6.3	63
Garlic	10	9.6	96
Ginger	10	9.1	91
Marjoram	10	< LOQ	< 25
Mustard seeds, black-brown	10	9.1	91
Mustard seeds, yellow	10	6.5	65
Nutmeg	10	10.5	105
Paprika powder, hot	10	7.1	71
Pepper, black	10	5.5	55

For some spices in table 10, recoveries below 60 % were measured due to the high polyphenol content and the mentioned disturbing effects. However, spices are usually only a minor component in final foods. Therefore, the spices were tested as samples containing 90 % cookie and 10 % of the respective spice. These samples were spiked again and recoveries are shown in the following table.

Table 11: Spices tested at 10 % level.

Sample composition	Expected concentration (mg/kg)	Concentration found (mg/kg)	Recovery (%)
Cookie + 10 % basil	10	7.4	74
Cookie + 10 % cloves	10	7.8	78
Cookie + 10 % marjoram	10	4.9	49
Cookie + 10 % pepper	10	7.3	73

Also cocoa was tested below 60 % in table 10. Since cocoa is mainly found in chocolate, several commercial chocolates with varying cocoa content were spiked and tested for recovery.

Table 12: Chocolates spiked with peanut.

Sample composition	Cocoa content (%)	Expected concentration (mg/kg)	Concentration found (mg/kg)	Recovery (%)
Chocolate 1	35	10	9.4	94
Chocolate 2	50	10	9.4	94
Chocolate 3	70	10	9.6	96

RIDASCREEN® FAST Peanut

Art. No. R6202

Samples from the market

Various samples have been collected from the market and analyzed. Presence of the allergen referring to the manufacturer's declaration was obtained for some samples. Please note that high sample dilutions lead to an increased measurement uncertainty. Furthermore, different processing conditions can have an impact on recovery.

Table 13: Measurement of samples from the market.

Sample description	Declaration by the manufacturer	Optical density	Concentration found (mg/kg)	Recovery (%)
Std 1		0.209		
Std 2		0.675		
Samples containing peanut (additional dilution factor of the extract)				
Cereal bar 1 (1:31000)	31 % peanuts	2.287	448900	145
Cereals (1:7000)	7 % peanuts	2.149	92384	132
Chocolate cream (1:25000)	25 % peanuts	1.198	109710	44
Cookie 1 (1:12000)	12 % peanuts	1.483	80106	67
Cookie 2 (1:21000)	21 % peanuts	1.855	221634	106
Cookie 3 (1:19000)	19 % peanuts	1.371	105644	56
Samples with possible traces of peanuts or nuts in general				
Cereals 2	May contain traces of peanut	0.683	2.5	
Cookie 4	May contain traces of peanut	1.112	4.0	
Cookie 5	May contain traces of peanut	0.236	< LOQ	
Cookie 6 (1:5)	May contain traces of peanut	4.085	> 100	
Chocolate 1	May contain traces of peanut	0.232	< LOQ	
Cereal bar 2	May contain traces of nuts	1.386	5.7	
Chocolate 1	May contain traces of peanut	0.234	< LOQ	
Ice cream 1	May contain traces of peanut	0.212	< LOQ	
Muesli 1	May contain traces of peanut	1.120	4.0	
Samples with no declaration of possibly contained peanuts				
Chocolate 2		0.238	< LOQ	
Chocolate 3		4.206	> 20	
Chocolate 4		0.184	< LOQ	
Cookies 7		0.446	< LOQ	
Ice cream 2		0.199	< LOQ	
Ice cream 3		0.180	< LOQ	
Muesli 2		0.217	< LOQ	
Rice crispies		0.183	< LOQ	



RIDASCREEN® FAST Peanut

Art. No. R6202

Additional information

Change of the concentrate of the Allergen extraction buffer (October 2012)


The old kit version contained a 20fold concentrate of the Allergen extraction buffer (AEB). Since October 2012 a new version contains a 10fold concentrate. Both buffers have to be diluted to the same final concentration of the buffer (1fold). The buffer composition was not changed and as, after the correct dilution, the buffer is the same in both cases, the same results were expected.

Before changing the format of the buffer concentrate a re-validation has been made under the supervision of AOAC Research Institute according to an agreed protocol. Negative matrices and spiked samples were measured to proof the assumption. Peanut free milk chocolate and breakfast cereals were extracted 30 times and all extracts were tested well below the LOQ (data not shown). Additionally, the matrices were also spiked with 5 mg/kg NIST peanut butter and extracted 30 times. All samples were detected positive.

Table 14: Results of two matrices (chocolate, cereals) spiked with a target value of 5 mg/kg peanut using the 10fold AEB for 30 extractions

	Milk chocolate spiked with 5 mg/kg peanut	Breakfast cereals spiked with 5 mg/kg peanut
Mean value (mg/kg peanut)	5.4	4.7
Recovery (%)	108	94
Standard deviation	0.37	0.37
CV (%)	6.9	7.9

Ruggedness of the test was investigated with small modifications to the regular test performance with regard to temperature (18 - 30 °C), incubation time, and volumes. The modifications applied did not significantly affect the test results (data not shown).



RIDASCREEN® FAST Peanut

Art. No. R6202

Results of collaborative studies

In addition, RIDASCREEN® FAST Peanut was validated with cereals, cookies, ice cream and chocolate in an independent study by three independent US laboratories under the supervision of AOAC Research Institute. All samples were prepared and distributed by FDA/CFSAN (College Park, MD). The FDA laboratory prepared 60 blank samples and 60 samples spiked at a level of 5 mg/kg peanut of each matrix and randomized the samples. Randomized sets of 20 blank and 20 spiked samples of each matrix were sent to each of the three independent laboratories. Samples have been prepared and analyzed according to a protocol of the AOAC Research Institute.

In table 15 a summary of the test results generated by the three independent laboratories is given. All of the spiked samples (240/240) were found positive (sensitivity = 100 %). Two out of 240 blank samples were measured positive (> 2.5 mg/kg peanut) corresponding to a specificity of 99.2 % (238/240). The mean value of all 240 spiked samples was 5.2 mg/kg with a recovery of 104 %.

Table 15: Summary of the test results analyzed by three independent laboratories

Target value	Found		% Agreement
	Positive	Negative	
Positive	240	0	100
Negative	2	238	99.2



RIDASCREEN® FAST Peanut

Art. No. R6202

Stability of the test

The stability of the test is routinely checked by R-Biopharm's quality assurance laboratory after defined storage intervals. Test kits are stored in a cold room at temperatures of 2 - 8 °C (35 - 46 °F). Before testing, the kit components are equilibrated to room temperature (20 - 25 °C / 68 - 77 °F).

Real time stability of the test is regularly controlled according to the total quality management schedule of the company.

Conclusion

With the described ELISA a sensitive and reliable ELISA is available which allows a quantitative and fast determination of the allergen or parts of it in food including all tested matrices.

Limitations


Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. Cross reactivities of the used antibodies have been determined for the pure food (e.g. corn flour). In a composed/processed food (e.g. corn bread) cross reactivities might be different. Interfering substances (e.g. polyphenols) can be detected by spike experiments.

For spiking, exemplary matrices/foods have been used.

The protein content and the protein composition may vary considerably between different peanut species. Therefore, different varieties may produce different results, since exemplary varieties were used for calibration.



RIDASCREEN® FAST Peanut


Art. No. R6202

Appendix 1 – QA Certificate

Quality Assurance Certificate

RIDASCREEN® FAST Peanut

REF Art. No.	R6202
LOT Lot	14515
Expiry	2017-05



	Lot. No.	Expiry		Lot. No.	Expiry
Microtiter plate	15425	2017-06	Red Chromogen Pro	14234	2018-04
Allergen extraction buffer	13435	2018-09	Stop solution	11255	2020-05
Standards	12505	2017-05			
Wash buffer	14495	2017-11			
Conjugate	12505	2017-05			

IFU (and other accompanying documents where applicable)

R-Biopharm AG, Darmstadt, Germany, certifies that this batch has been approved by the Quality Assurance Department and conforms to specifications.

Please note: The absorbance for the standards may decrease during the shelf life of the kit. The general shape of the curve will remain similar, while the slope might change slightly. Furthermore refer to product leaflet section 8. Indication of instability or deterioration of reagents.

R-Biopharm AG, Quality Assurance Department

Date: 2015-12-17



Remark: This document is created electronically and therefore valid without a signature.

R-Biopharm AG
 Postal address: An der neuen Bergstraße 17
 64297 Darmstadt, Germany
 Corporate seat: Pfungstadt
 Phone: +49 (0) 61 51 - 81 02-0
 Fax: +49 (0) 61 51 - 81 02-40
 www.r-biopharm.com

Chairman of supervisory board: Dietrich Müllat
Board of directors:
 Dr. Ralf M. Döhler (CEO),
 Dr. Carsten Bruns,
 Jochen Hirsch,
 Dr. Peter Schubert

Reg.No.: Amtsgericht Darmstadt, HRB 8321
Tax No.: 07242 / 02-888
VAT Reg. No.: DE 111 657 409





RIDASCREEN® FAST Peanut

Art. No. R6202

Appendix 2 – Spike instructions

1. Spike with NIST RM® 2387 peanut butter

The National Institute of Standards and Technology (NIST) in the USA offers a Standard Reference Material® 2387 Peanut butter which has a concentration of 90 % peanut.

- Dissolve 22 mg of NIST peanut butter in 20 ml Allergen extraction buffer and shake for 2 hours upside down or by a rotator (concentration: 1 mg/ml peanut).
- Dilute 1:10 with Allergen Extraction buffer (e.g. 500 µl of stock solution + 4500 ml of Allergen extraction buffer) and vortex thoroughly (concentration: 100 µg/ml peanut).
- Add 100 µl of this solution to 1 g of sample to obtain a 10 mg/kg spike.
- Extract the sample as described in the instructions for use of RIDASCREEN® FAST Peanut.

2. Spike with American Peanut Butter

American Peanut Butter from the supermarket (e.g. Barney's Best Peanut Butter) contains approx. 90 % of roasted peanut.

- Dissolve 22 mg of peanut butter in 20 ml Allergen Extraction buffer and shake for 2 hours upside down or by a rotator (concentration: 1 mg/ml peanut).
- Dilute 1:10 with Allergen extraction buffer (e.g. 500 µl of stock solution + 4500 ml of Allergen Extraction buffer) and vortex thoroughly (concentration: 100 µg/ml peanut).
- Add 100 µl of this solution to 1 g of sample to obtain a 10 mg/kg spike.
- Extract the sample as described in the instructions for use of RIDASCREEN® FAST Peanut.



RIDASCREEN® FAST Peanut

Art. No. R6202

Appendix 3 – Equipment

Pipettes:

- 10 - 100 µl Eppendorf Research Plus, Eppendorf
- 20 - 200 µl Eppendorf Research Plus, Eppendorf
- 100 - 1000 µl Eppendorf Research Plus, Eppendorf
- Multichannel: 30 - 300 µl Eppendorf Research Plus, Eppendorf

Tips:

- epT.I.P.S.® Standard/Bulk 2 - 200 µl, Eppendorf
- epT.I.P.S.® Standard 50 - 1000 µl, Eppendorf

Multistepper:

- Multipette® Xstream Eppendorf Research Plus, Eppendorf

Multistepper tips:

- 2.5; 5 and 50 ml Combitips Advanced®, Eppendorf

Serological pipette:

- 5 ml and 25 ml CELLSTAR® Serological Pipette, Greiner Bio-One

Pre-plate:

- Mikrotiter Assembly Breakable Strip 1x8, Thermo Scientific
- Low binding from Greiner Bio-One, Cat.-No. 655101

Photometer:

- Tecan Sunrise, Tecan

Wash:

- Dispensette®, Brand
- Brand Tech® 8-channel Manifold, Brand

For further information or orders please contact R-Biopharm AG:

International Sales:

Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-40
E-mail: sales@r-biopharm.de

Order Department:

Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-mail: orders@r-biopharm.de