RIDA[®]QUICK Soya Extraction buffer

Art. No. R7113

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RIDA[®]QUICK Soya Extraction buffer

1. Intended Use

The RIDA[®]QUICK Soya Extraction buffer is used for the sample preparation for the determination with RIDA[®]QUICK Soya (Art. No. R7103).

Related products

RIDA[®]QUICK Soya (R7103)

Set of 3 processed Soya Assay Controls (R7132)

2. Content

Art. No.: R7113 (2 x 100 ml)

3. Materials required but not provided

- centrifuge + centrifugal vials
- shaker
- water bath
- laboratory mincer/grinder, mortar, ultra-turrax or homogenizer
- fluted filter
- graduated pipettes

4. Warnings and precautions for the users

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

This kit may contain hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

5. Storage instructions

Store the kit at room temperature. **Do not freeze any test kit components.**

6. Preparation of Samples

Working devices such as a mill, glass vials or spatulas must be cleaned before and after each sample preparation to remove any remains of allergen and to avoid contamination.

- Grind at least 5 g of the sample carefully and mix thoroughly
- weigh 1 g of the sample and add 7.5 ml RIDA[®]QUICK Soya Extraction buffer (the extraction buffer should already have a temperature of approx. 60 °C (140 °F))
- mix intensively and incubate for 10 min at 100 °C (212 °F) in a boiling water bath, afterwards cool down
- centrifuge: 10 min, at least 2500 g, if possible at 4 °C (39 °F) and/or filter the extract (alternatively 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- continue as described in the instruction for use of RIDA[®]QUICK Soya (Art. No. R7103) point 8. Test implementation (swab test), subitem 8

In general

Samples tested negative still could contain an allergen contamination below the limit of detection of the assay, or they might contain other allergen components like lipids for example.

Due to the multitude of food types, matrix effects cannot be excluded. In processed food (e.g. heat treatment, dehydration, etc.), proteins may be altered or fragmented, this may have an impact on the recovery/cross reactivity.

For evaluation of the cross reactivity only one exemplary sample was analyzed, other samples may show a different result. All cross reactivities and exemplary analyzed matrices are described in the Validation report

Recommendation

In order to ensure a high analytical performance we recommend:

- to use also allergen-free and allergen- containing (spiked) samples as test controls
- due to the multitude of food types, matrix effects cannot be excluded. To ensure an accurate result spike experiments are recommended
- for quantification of the results use the ELISA RIDASCREEN[®]FAST Soya (Art. No. R7102)
- to perform SureFood[®] PCR for confirmation of the result

For further product information and application notes, please contact sales@r-biopharm.de.

The product information folder with further information is available on request from your local distributor or R-Biopharm AG.

Further application notes

–Extraction of the Set of 3 processed Assay Controls (Art. No. R7132)
for the RIDA®QUICK Soya

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG