RIDASCREEN[®] FAST Sesame

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RIDASCREEN[®] FAST Sesame

Brief information

RIDASCREEN[®]FAST Sesame (Art. No. R7202) is a sandwich enzyme immunoassay for the quantitative analysis of sesame or parts of sesame in food like cereals, pastries, ice cream and chocolate.

All reagents required for the enzyme immunoassay - including standards - are contained in the test kit.

The test kit is sufficient for 48 determinations (including standards).

A microtiter plate spectrophotometer is required for quantification.

Sample preparation:	homogenization, extraction and centrifugation
Time requirement:	sample preparation (for 10 samples)approx. 20 min test implementation (incubation time)30 min
Limit of detection:	0.24 mg/kg (ppm) sesame corresponding to 0.000024 % sesame
Limit of quantification:	2.5 mg/kg sesame corresponding to 0.00025 % sesame
Specificity:	The antibody specifically detects proteins from sesame. No cross-reaction with cashew nut, brazil nut, pecan, hazelnut, coconut, macadamia nut, walnut, chest nut, sunflower seeds and lima beans.

The specificity of the RIDASCREEN[®]FAST Sesame test was determined by analyzing the cross-reactivities to corresponding substances in buffer system. In samples, the specificity may deviate from those determined in the buffer system due to matrix effects. Prior to the analysis of cross-reactive substances, the user has to determine the Limit of Detection and the Recovery for the substance in the respective sample matrix. The test cannot discriminate between analytes and cross-reactive substances.

Related products:

bioavid Lateral Flow Sesame (Art. No. BL609-10/-25)

SureFood® ALLERGEN ID Sesame (Art. No. S3108)

SureFood® ALLERGEN QUANT Sesame (Art. No. S3208)

1. Intended use

RIDASCREEN®FAST Sesame (Art. No. R7202) is a sandwich enzyme immunoassay for the quantitative analysis of sesame or parts of sesame (as ingredient components or contamination) in raw or heated food like cereals, pastries, ice cream and chocolate.

2. General

Sesame seeds (*Sesamum indicum*) are used as raw or roasted ingredients in pastries, ice cream, cereals, chocolate or for the production of oil or paste (Tahina). Labeling of sesame is mandatory in the European Union, Canada and Australia.

Sesame seed food allergy is an increasingly recognized health burden, the symptoms vary from allergic eczema to life threatening anaphylactic reactions. Until now, seven different proteins of sesame were identified as allergens (Ses i 1-7). The only effective treatment to protect sensitized patients from allergic symptoms is a strict avoidance of sesame containing food.

3. Test principle

The basis of the test is the antigen-antibody reaction. The wells of the microtiter strips are coated with specific antibodies to sesame-protein. By adding the standard or sample solution to the wells, present sesame-protein will bind to the specific capture antibodies. The result is an antibody-antigen-complex. Components not bound by the antibodies are then removed in a washing step. Then antibody conjugated to peroxidase is added. This antibody-conjugate is bound to the Ab-Ag-complex. An antibody-antigen-antibody (sandwich) complex is formed. Any unbound enzyme conjugate is then removed in a washing step. Enzyme substrate and chromogen are added to the wells and incubated. Bound enzyme conjugate converts the chromogen into a blue product. The addition of the stop reagent leads to a color change from blue to yellow. The measurement is made photometrically at 450 nm. The absorbance is proportional to the sesame concentration of the sample.

4. Reagents provided

Each kit contains sufficient materials for 48 measurements (including standard analyses). Each test kit contains:

- 1 x Microtiter plate with 48 wells (6 strips with 8 removable wells each)
coated with anti-sesame antibodies
- 5 x Standards ^{*)}, 1.3 ml each
0 ppm (zero standard), 2.5 ppm, 5 ppm, 10 ppm, 20 ppm sesame
in aqueous solution, ready to use
- 1 x Conjugate (7 ml).....red cap
peroxidase conjugated antibody, ready-to-use
- 1 x Red Chromogen Pro (10 ml) brown cap
substrate/chromogen solution, stained red
- 1 x Stop solution (14 ml) yellow cap
contains 1 N sulfuric acid
- 1 x Allergen Extraction buffer (125 ml)green cap
as a 10fold concentrate
- 1 x Washing buffer (100 ml) brown cap
as a 10fold concentrate

*) The dilution factor 20 for the sample has already been considered when labeling. Therefore, the sesame concentrations of samples can directly be read from the standard curve.

5. Reagents required but not provided

5.1. Equipment:

- microtiter plate spectrophotometer (450 nm)
- centrifuge + centrifugal glass vials
- water bath
- laboratory mincer / grinder, pestle and mortar, Ultra-Turrax or homogenisator
- graduated pipettes
- paper filter
- variable 20 µl - 200 µl and 200 - 1000 µl micropipettes
- if possible multichannel pipette and/or multistepper pipette

5.2. Reagents:

- distilled or deionized water
- skim milk powder (SMP, food quality)

6. Warnings and precautions

This test should only be carried out by trained laboratory employees. The instruction for use must be strictly followed.

The washing buffer contains thimerosal. Avoid skin contact.

The stop solution contains 1 N sulfuric acid (R36/38, S2-26).

This kit may contain further hazardous substances. For hazard notes on the contained substances please refer to the appropriate material safety data sheets (MSDS) for this product, available online at www.r-biopharm.com.

7. Storage instructions

Store the kit at 2 - 8 °C (35 - 46 °F). Do not freeze any test kit components.

Return any unused microwells to their original foil bag, reseal them together with the desiccant provided and further store at 2 - 8 °C (35 - 46 °F).

Red Chromogen Pro is light sensitive, therefore, avoid exposure to direct light.

No quality guarantee is accepted after the expiry date on the kit label.

Do not interchange individual reagents between kits of different lot numbers.

8. Indication of instability or deterioration of reagents

- any bluish coloration of the red stained Red Chromogen Pro prior to test implementation
- a value of less than 0.8 absorbance units ($A_{450\text{ nm}} < 0.8$) for standard 5

9. Preparation of Samples

The **Allergen Extraction buffer** is provided as a **10fold concentrate**. Before dilution of the buffer concentrate dissolve any crystals in a water bath at 37 °C (99 °F) completely and mix well. After that dilute the heated buffer concentrate 1:10 (1+9) with distilled water before use (e.g. 100 ml buffer concentrate + 900 ml distilled water, or alternatively use the whole content of the bottle and fill up to

1250 ml). Add 5 g skim milk powder (SMP) per 100 ml diluted Allergen Extraction buffer (respectively 62.5 g to 1250 ml). The 1:10 diluted Allergen Extraction buffer without SMP is stable for approx. twelve weeks and with SMP for one day at 2 - 8 °C (35 - 46 °F).

9.1. Test preparation

Unclean laboratory equipment can result in a contamination of samples or test components with sesame. Therefore, gloves should be worn during the test procedure and the following recommendations should be applied:

- thoroughly clean surfaces, glass vials, mixers and additional equipment before and after each sample
- sample extraction and ELISA-procedure should be performed in separated rooms
- thoroughly pound a representative amount of the sample (5 - 50 g), crumble and homogenize it

9.2. Solid samples

- weigh 1 g of a representative sample and add 20 ml preheated (60 °C / 140 °F) and finally diluted Allergen Extraction buffer with SMP, close the vial, mix vigorously and incubate the solution for 10 min at 60 °C / 140 °F in a water bath
- let the sample cool down (e.g. in an ice bath) and centrifuge 10 min 2500 g if possible at 4 °C (39 °F) and / or filter the extract (alternatively, 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl extract per well in the assay

9.3. Liquid samples

- take 1 ml of a representative sample and add 19 ml preheated (60 °C / 140 °F) and finally diluted Allergen Extraction buffer with SMP, close the vial, mix vigorously and incubate the solution for 10 min at 60 °C / 140 °F in a water bath
- let the sample cool down (e.g. in an ice bath) and centrifuge 10 min 2500 g if possible at 4 °C (39 °F) and / or filter the extract (alternatively, 2 ml of the extract can be centrifuged at high speed for 10 min in reaction caps by using a microcentrifuge)
- use 100 µl extract per well in the assay

Remark:

Further dilutions should be prepared with the finally diluted Allergen Extraction buffer with SMP.

The sample extracts can be stored at 2 - 8 °C (35 - 46 °F) for 3 days. Unused extracts can be stored at -20 °C (-4 °F) for several months.

10. Test implementation

10.1. Preliminary comments

Bring all reagents to room temperature (20 - 25 °C / 68 - 77 °F) before use.

The **Washing Buffer** is provided as a 10fold concentrate. Before use, the buffer has to be diluted 1:10 (1+9) with distilled water (i.e. 100 ml buffer concentrate + 900 ml dist. water). Prior to dilution, dissolve eventually formed crystals by incubating the buffer in a water bath at 37 °C (99 °F). The diluted buffer is stable at 2 - 8 °C (35 - 46 °F) for approx. four weeks.

10.2. Test procedure

Carefully follow the recommended washing procedure. Do not allow microwells to dry between working steps.

Do not use more than three strips (24 wells) at the same time. In the case of more than three strips, a second uncoated plate (e.g. low binding from Greiner bio-one, Cat.-No. 655101) should be used to avoid a time shift over the microtiter plate. All standards and samples (at least 150 µl for every well) are pipetted into the uncoated plate and then 100 µl are quickly transferred to the coated microtiter plate with an 8-channel pipette.

1. Insert a sufficient number of wells into the microwell holder for all standard and samples to be run in duplicate. Record standard and sample positions.
2. Add 100 µl of each standard solution or prepared sample to separate duplicate wells and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
3. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1). Repeat two more times.

4. Add 100 µl of the enzyme conjugate to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F).
5. Pour the liquid out of the wells and tap the microwell holder upside down vigorously (three times in a row) against absorbent paper to ensure complete removal of liquid from the wells. Fill all the wells with 250 µl washing buffer (see 10.1.) and pour out the liquid again. Repeat two more times.
6. Add 100 µl Red Chromogen Pro to each well. Mix gently by shaking the plate manually and incubate for 10 min at room temperature (20 - 25 °C / 68 - 77 °F) in the dark.
7. Add 100 µl of the stop solution to each well. Mix gently by shaking the plate manually and measure the absorbance at 450 nm. Read within 10 minutes after addition of stop solution.

11. Results

A special software, the RIDA[®]SOFT Win / RIDA[®]SOFT Win.net (Art. Nr. Z9996), is available for evaluation of the RIDASCREEN[®] enzyme immunoassays. The calculation should be done by use of a cubic spline function. The course of the standard curve is shown in the Quality Assurance Certificate enclosed in the test kit.

In comparison with the certificate, higher values of the absorbance ($A_{450\text{ nm}}$) for the calibration curve, especially for the zero standard, may be a result of insufficient washing or sesame contamination.

A further dilution and new detection of the samples is recommended for absorbance values ($A_{450\text{ nm}}$) > standard 5.

Please note:

When working in accordance with the regulation stated, the dilution factor is 20. The sesame concentration can be read directly from the standard curve (see 4. *)
- the sample dilution factor of 20 is already taken into account.

For sample dilutions of more than 1:20, the further dilution factor must be considered for the calculation of the sesame concentration.

In general:

Due to the multitude of food types, matrix effects cannot be excluded.

In processed food, proteins may be altered or fragmented, this may have an impact on the recovery.

Remark:

Samples tested negative could still contain a sesame contamination below the limit of detection of the assay, or they might contain other protein-free sesame components (e.g. oil).

Recommendation:

In order to ensure a high analytical performance:

- each sample material should be analyzed in duplicates
- employ also sesame protein free and sesame protein containing (spiked) samples as test controls

The data corresponds to our present state of technology and provides information on our products and their uses. R-Biopharm makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. Defective products will be replaced. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. R-Biopharm shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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are registered trademarks of R-Biopharm AG