**SureFood® ALLERGEN QUANT Sesame
(100 Reakt.)**

Art.-No. S3208

Version 1.3

Example:

$C_{t \text{ sample}}$:	28.5	$X \text{ [ppm]} = \frac{10 \cdot \left(\frac{28.5 - 40.54}{-3.251} \right)}{\left(\frac{29.2 - 40.54}{-3.251} \right)} * 40 = 65.7 \text{ ppm}$
$C_{t \text{ QUANTARD}}$:	Ø 29.2	
a:	40.54	
s:	- 3.251	
c:	40	

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Version 1.3

4. Setup

	Blockcycler/LC 480	LightCycler® / Rotor-Gene Q
Initial Denaturation (HOLD) Cycles	5 min, 95°C 45	1 min, 95°C 45
Denaturation Annealing/Extension (CYCLE)	15 sec, 95°C 30 sec, 60°C	10 sec, 95°C 15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Reporter Dye: FAM Quencher Dye: BHQ Passive reference: none	LightCycler® Channel: 530 or F1 Acquisition mode: Single in extension phase Rotor-Gene Q Reporter Dye: Green

Interpretation of results

The negative and positive controls need to give the correct results.

A sample is stated **negative**, if the sample DNA shows no amplification in the detection system or the calculated result is below the limit of detection. In case of a **negative** result the inhibition control of the sample must be **positive**. If this is not the case, the sample contains PCR-inhibiting substances. Under these circumstances an evaluation of the samples is not possible. DNA isolation and purification for the sample need to be improved.

For the calculation of sample DNA with an amplification curve mark the standards, the controls and the samples and make the evaluation according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The value for the slope of the standard curve has to be between -3.2 and -3.6 and the correlation coefficient > 0.98. If the standard curve delivers different values it should not be used for an evaluation.

For a final calculation of mg allergenic ingredient/kg food sample following formula is introduced:

$$X [ppm] = \frac{10^{\left(\frac{Ct_{sample} - a}{s}\right)}}{10^{\left(\frac{Ct_{QUANTARD} - a}{s}\right)}} * c_{QUANTARD}$$

a = axis intercept
s = slope
c = concentration
Ct_{QUANTARD} = Ø of double assay

SureFood® ALLERGEN QUANT Sesame (100 React.)

Art.-No. S3208

Version 1.3

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Sesame Reaction Mix / Inhibition Control Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2. Preparation of the standard DNA dilutions

Dilute the Standard DNA (**Code 5**) in 1:10 steps in Dilution Buffer (**Code 4**) in order to prepare different concentrations. Prepare a dilution series of 4 steps. Prepare 4 reaction tubes (labeled S1 to S4) and add 45 µl Dilution Buffer each. The following procedure is recommended:

Standard	dilutions	copy number per µl	final copy number per reaction*
S1	45 µl Dilution Buffer + 5 µl Standard DNA	10.000 copies	50.000 copies
S2	45 µl Dilution Buffer + 5 µl DNA of S1	1000 copies	5.000 copies
S3	45 µl Dilution Buffer + 5 µl DNA of S2	100 copies	500 copies
S4	45 µl Dilution Buffer + 5 µl DNA of S3	10 copies	50 copies

***Note:** 5 µl of standard DNA are used for each calibration point. The final copy number per reaction is to be entered in the analysis software of the real-time PCR detection system.

If the diluted DNA standards (S1 to S4) are not immediately used, store them at -20°C. The dilution series are stable up to two months at -20°C. Before use allow the reagents to thaw, mix them on a vortex and centrifuge carefully before opening and use.

3. Preparation of the PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells or capillaries.
- Close the tube of the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes.
- Pipette 5 µl of the SureFood® QUANTARD Allergen 40 DNA and the standard dilutions into the designated tubes.
- Centrifuge all tubes/plates or capillaries shortly at low speed.
- Place tubes/plates into the PCR instrument and start the run according to the setup.

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Version 1.3

Description

The test detects the quantitative DNA amount of sesame in food samples according to EU regulation 2007/68 EC. For the quantitative determination the use of the laboratory reference material SureFood® QUANTARD Allergen 40 containing 40 mg white sesame seeds/ kg food sample and a standard dilution is recommended. The real-time PCR assay can be used with established real-time PCR instruments.

Remarks

In general this procedure is applicable for all homogenous samples. The quantitative determination is inapplicable for swabs, sponges and inhomogeneous fluids/liquids. The method is validated for quantitative determination of 0.5 to 400 mg allergenic ingredient / kg food sample using the reference material SureFood® QUANTARD Allergen 40 containing 40 ppm allergenic ingredient / kg food sample.

Limit of Detection

The SureFood® ALLERGEN QUANT Sesame PCR has a limit of detection of ≤ 0.4 ppm. The limit of quantification is 0.5 ppm using SureFood® PREP Allergen. The limit of detection and quantification were determined using the SureFood® QUANTARD Allergen 40 (matrix: corn flour). The assay limit of detection depends on sample matrix, processing grade, DNA-preparation and DNA-content.

DNA-preparation

For DNA-preparation the use of SureFood® PREP Allergen is recommended.

Kit components and storage

2x	Sesame Reaction Mix (1.0 ml)	(Code 1)
2x	Inhibition Control Mix (1.0 ml)	(Code 2)
1x	Taq Polymerase (20 µl)	(Code 3)
1x	Dilution Buffer (2.0 ml)	(Code 4)
1x	Standard DNA (80 µl)	(Code 5)

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- **SureFood® QUANTARD Allergen 40 Art.-No.: S3301**
- Real-time PCR instrument
- Real-time PCR consumables (plates, tubes, foils, capillaries)
- Pipettes with filter tips

Protocol

1. Preparation of the master-mix

Calculate the total number of reactions needed (samples, controls and standards).

Reactions needed for the sesame detection:

For each run: 4 reactions for the standard curve
4 reactions for controls (1x no-template control, 1x negative extraction control and 2x SureFood® QUANTARD Allergen 40 DNA)

For each sample: at least 1 reaction for each sample DNA