**SureFood® ALLERGEN 4plex Soya/Celery/Mustard+IAC
(100 Reakt.)**

Art. Nr. S3401

Version 1.4

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Description

The multiplex test detects DNA of soya, celery and mustard according to directive (EC) 1169/2011. The real-time PCR assay can be used with established real-time PCR instruments, equipped for the detection of four fluorescence emissions at 510 nm, 580 nm, 610 nm and 660 nm (FAM, VIC/HEX, ROX and Cy5) at the same time. The technical validation was performed on Agilent Mx3005P, Roche LightCycler® 480 II*, Applied Biosystems 7500, BioRad CFX96 und Qiagen Rotor-Gene Q.

The kit also contains an internal amplification control (inhibition control, IAC) to examine DNA for inhibiting substances which interfere with the PCR mix. The Reaction Mix contains all reagents necessary for PCR as well as an artificial target (concentration: 500 copies / reaction). The target is amplified during the PCR run and detected in the VIC/HEX-channel. When the DNA contains PCR inhibiting substances, amplification of the target will be suppressed and the signal will be affected.

Some examples for PCR inhibiting substances are alcohols (e.g. ethanol, isopropanol), surfactants (e.g. CTAB, SDS, Triton X100) and salts (e.g. sodium chloride). Also spices, herbs, algae, cocoa and further sample matrices might have PCR inhibiting effects.

Limit of Detection

The SureFood® ALLERGEN 4plex Soya/Celery/Mustard+IAC PCR has a limit of detection of ≤ 0.4 mg/kg allergenic substance in non-processed corn flour using SureFood® PREP Advanced, protocol 1.

The limit of detection of the complete assay (DNA extraction and real-time PCR) depends on sample matrix, processing grade, DNA-preparation and DNA-content.

Very high levels of soya or mustard might influence the sensitivity of the celery detection system.

DNA-preparation

For DNA-preparation the use of SureFood® PREP Advanced, protocol 1 is recommended.

Kit components and storage

2x	Reaction Mix (1.1 ml)	(Code 1)
1x	Taq Polymerase (11 µl)	(Code 2)
1x	Positive Control (200 µl)	(Code 3)

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- real-time PCR instrument, equipped with four detection channels (510 nm, 580 nm, 610 nm, 660 nm)
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips
- unpowdered disposable gloves
- Vortex mixer
- micro centrifuge with a rotor for the reaction tubes

* Note: For users of the Roche LightCycler® 480 II instrument a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art. No. F4009) must be used for the color compensation of such devices.

Protocol

1. Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, extraction control, positive control. Beside the soya/celery/mustard detection assays the Reaction Mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use. The tube of the Taq Polymerase should be kept at -20°C and not be mixed by vortexing.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10% excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix each master-mix well and centrifuge shortly before use.

2. Setup

	Blockcyler*	Rotorcyler
Initial Denaturation (HOLD) Cycles	5 min, 95°C 35	1 min, 95°C 35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Detection of celery: in the FAM-channel (LC480 II channel 465-510) Detection of the internal amplification control: in the VIC/HEX-channel (LC480 II channel 533-580) Detection of mustard: in the ROX-channel (LC480 II channel 533-610) Detection of soya: in the Cy5-channel (LC480 II channel 618-660)	Detection: End of extension phase Detection of celery: in the FAM-channel (Green) Detection of the internal amplification control: in the VIC/HEX-channel (Yellow) Detection of mustard: in the ROX-channel (Orange) Detection of soya: in the Cy5-channel (Red)
Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage: http://www.congen.de/en/company/downloads		

* Note: For users of the Roche LightCycler® 480 II instrument a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art. No. F4009) must be used for the color compensation of such devices.

3. Preparation of the real-time PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes/wells.
- Close the negative control (the negative control is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes/wells and close them.
- Pipette 5 µl of Positive Control into the designated tubes/wells and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the real-time PCR instrument and start the run according to the setup.

Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer.

The control reactions need to give the correct results.

Soya DNA is detected in the Cy5-channel, mustard DNA is detected in the ROX-channel, celery DNA is detected in the FAM-channel and the internal amplification control is detected in the VIC/HEX-channel.

A sample is stated **positive** for the respective parameter (soya/celery/mustard), if the sample DNA shows amplification in the respective channel (Cy5/FAM/ROX).

A sample is stated **negative** for the respective parameter (soya/celery/mustard), if the sample DNA shows no amplification in the respective channel (Cy5/FAM/ROX) and the internal amplification control (VIC/HEX) of the sample is **positive** (see also table on page 7).

The internal amplification control is a heterologous PCR system. It can be influenced by other factors than the specific detection system.

result in the respective channel				
FAM-channel celery	ROX-channel mustard	Cy5-channel soya	VIC/HEX-channel amplification control	result
positive	negative	negative	positive/ negative	celery DNA detected
negative	positive	negative	positive/ negative	mustard DNA detected
negative	negative	positive	positive/ negative	soya DNA detected
negative	negative	negative	negative	invalid [#]

[#]In case of a negative result in all channels including the VIC/HEX-channel, the sample contains PCR inhibiting substances. Under these circumstances an evaluation of the sample is not possible. DNA isolation and purification for the sample need to be improved.

Product Information

- Validation Report

Technical Support

For further questions please contact your distributor or send an e-mail to sales@r-biopharm.de.

Distribution and ordering

R-Biopharm AG
An der neuen Bergstrasse 17,
64297 Darmstadt, Germany
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-Mail: orders@r-biopharm.de
www.r-biopharm.com