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Interpretation of results

The evaluation has to be made according to the usual analysis program recommended by the real-time PCR instrument manufacturer. The control reactions need to give the correct results.

A sample is stated **positive** for the respective parameter (peanut/hazelnut/walnut), if the sample DNA shows amplification in the respective channel. A sample is stated **negative** for the respective parameter (peanut/hazelnut/walnut), if the sample DNA shows no amplification in the respective channel and the internal signal (VIC/HEX-channel) of the sample is **positive** (see also table on page 6).

result in the respective channel				result
FAM-channel peanut	ROX-channel walnut	Cy5-cannel hazelnut	VIC/HEX-channel amplification control	
positive	negative	negative	positive/ negative	peanut DNA detected
negative	positive	negative	positive/ negative	walnut DNA detected
negative	negative	positive	positive/ negative	hazelnut DNA detected
negative	negative	negative	negative	invalid [#]

[#] In case of a negative result in all channels including the VIC/HEX-channel, the sample contains PCR inhibiting substances. Under these circumstances an evaluation of the sample is not possible. DNA isolation and purification for the sample need to be improved.

Product Information

- Validation Report
- Product Information
- Material Safety Data Sheet

Technical Support

For further questions please send an e-mail to sales@r-biopharm.de.

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Description

This multiplex test detects specific DNA sequences of peanut, hazelnut and walnut. The real-time PCR assay can be used with established real-time PCR instruments, equipped for the detection of four fluorescence emissions at 510 nm, 580 nm, 610 nm and 660 nm (FAM, VIC/HEX, ROX and Cy5) at the same time (Agilent MxSeries, Roche LightCycler® 480 II†, ABI PRISM, BioRad CFX96 etc.).

Limit of Detection

The SureFood® ALLERGEN 4plex Peanut/Hazelnut/Walnut+IAC real-time PCR has a limit of detection of ≤ 1 mg/kg. The assay limit of detection depends on sample matrix, processing grade, DNA preparation and DNA content.

DNA-preparation

For DNA-preparation the use of SureFood® PREP Allergen is recommended.

Kit components and storage

2x	Reaction Mix (1.1 ml)	(Code 1)
1x	Taq Polymerase (11 µl)	(Code 2)
1x	Positive Control (200 µl)	(Code 3)

Store all reagents at -20°C and protected from light.

Additionally required equipment and materials

- real-time PCR instrument, equipped with four detection channels
- real-time PCR consumables (plates, tubes, foils, caps)
- pipettes with filter tips

Protocol

1. Preparation of the master-mix

Calculate the total number of reactions needed (samples and control reactions). Recommended control reactions: negative control, positive control and extraction control for each sample. Beside peanut/hazelnut/walnut detection assays the master-mix includes an internal amplification control (inhibition control) for each reaction.

It is also recommended to prepare the master-mix with 10 % additional volume in order to compensate reagent loss. Allow the reagents to thaw, mix by vortexing and centrifuge before opening and use.

Example for the calculation and preparation of 10 reactions:

Components for master-mix	Amount per reaction	10 reactions (with 10 % excess)
Reaction Mix	19.9 µl	218.9 µl
Taq Polymerase	0.1 µl	1.1 µl
Total volume	20 µl	220 µl

Mix the master-mix well and centrifuge shortly before use.

† Note: For users of the Roche LightCycler® 480 II instrument a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art.-No. S4009) is recommended for the color compensation of such devices.

2. Setup

	Blockcycler / LightCycler® 480 II†	Rotor-Gene Q
Initial Denaturation (HOLD) Cycles	5 min, 95°C 35	1 min, 95°C 35
Denaturation	15 sec, 95°C	10 sec, 95°C
Annealing/Extension (CYCLE)	30 sec, 60°C	15 sec, 60°C
Temperature Transition Rate/ Ramp Rate	Maximum	Maximum
Fluorescence Detection Setup	Detection: End of extension phase Detection of peanut: in the FAM-channel (LC480 II channel 465-510) Detection of the internal amplification control: in the VIC/HEX-channel (LC480 II channel 533-580) Detection of walnut: in the ROX-channel (LC480 II channel 533-610) Detection of hazelnut: in the Cy5-channel (LC480 II channel 618-660)	Detection: End of extension phase Detection of peanut: in the FAM-channel (Green) Detection of the internal amplification control: in the VIC/HEX-channel (Yellow) Detection of walnut: in the ROX-channel (Orange) Detection of hazelnut: in the Cy5-channel (Red)

Detailed information on the setup of several real-time PCR devices is available at the CONGEN homepage:

<http://www.congen.de/en/information/device-settings-for-surefood-real-time-pcr-kits.html>

3. Preparation of the PCR-mix

- Pipette 20 µl of the master-mix into appropriate tubes or wells.
- Close the tubes/wells of the negative control (it is ready for PCR without any addition).
- Pipette 5 µl of sample DNA into the designated tubes and close them.
- Pipette 5 µl of the Positive Control into the designated tubes and close them.
- Centrifuge all tubes/plates shortly at low speed.
- Place tubes/plates into the PCR instrument and start the run according to the setup.

† Note: For users of the Roche LightCycler® 480 II instrument a Color Compensation is necessary. The SureCC Color Compensation Kit I (Art.-No. S4009) is recommended for the color compensation of such devices.