






RIDA® QUICK Zearalenon ECO

Art. No. R5506

Please see the instructions for use (IFU) and the reader manuals for details

Bring all components to room temperature (20 - 25 °C/68 - 77 °F) and perform the entire testing process at room temperature. For evaluation, use the RIDA®SMART APP software (Art. No. ZRSAM, see IFU chapter 10) in combination with a validated Android smartphone or with the RIDA®SMART BOX (Art. No. ZRSA-SB) and a smartphone.







These brief instructions have been prepared with the greatest care, but are intended only as a guide and brief information and do not replace the official instructions for use (IFU) or the careful reading of these. Only the information in the instructions for use applies.

RIDA® QUICK Zearalenon ECO

Art. Nr. R5506

Details: siehe Durchführungsanweisung (TKB) und Auswertegerät/-software-Handbücher

Alle verwendeten Komponenten vor dem Einsatz im Test auf Raumtemperatur bringen (20 - 25 °C) und den gesamten Prozess bei Raumtemperatur durchführen. Die Auswertung erfolgt mit der RIDA®SMART APP Software (Art. Nr. ZRSAM, siehe Kapitel 10 der Anleitung) in Kombination mit einem validierten Android Smartphone oder mit der RIDA®SMART BOX (Art. Nr. ZRSA-SB) und einem Smartphone.



Diese Kurzanleitung wurde mit größter Sorgfalt erstellt, dient jedoch nur als Hilfestellung und Kurzinformation und ersetzt nicht die offizielle Durchführungsanweisung (TKB) sowie deren aufmerksames Lesen. Es gelten ausschließlich die Angaben in der Durchführungsanweisung.